

## MUELLER-HINTON II AGAR

An antimicrobial susceptibility testing medium which fits the requirements of EUCAST. The medium has extremely low concentrations of thymine and thy-midine as well as appropriate levels of calcium and magnesium ions.



Dehydrated media	
Code number:	500 g: MHT20500, 5 kg: MHT25000
Colour:	Yellowish
Appearance:	Homogeneous hygroscopic powder
pH before autoclaving (25 °C):	7,1 – 7,5

**Direction for Mueller-Hinton II Agar:** Suspend **38 g** in one litre of distilled water and heat with frequent agitation until the medium boils well. Sterilise by autoclaving at 121 °C for 15 minutes.

Direction for Mueller-Hinton II Blood Agar, EUCAST: Suspend 38 g in one litre of distilled water and heat with frequent agitation until the medium boils well. Sterilise by autoclaving at 121 °C for 15 minutes. Cool to 50 °C and add aseptically 50 ml of sterile defibrinated horse blood and 0,02 g  $\beta$ -NAD. Mix well before pouring.

**Direction for Mueller-Hinton II Chocolate Agar:** Suspend **19 g** in 460 ml of distilled water and heat with frequent agitation until the medium boils well. Sterilise by autoclaving at 121 °C for 15 minutes. Cool to 50 °C and add aseptically **35 ml of sterile defibrinated blood** and "chocolate" by heating at 80 °C for 10 min. Cool to 50 °C. Dissolve the contents of **one vial of Growth Factor Mixture Hydration Fluid** with 5 ml of sterile distilled water and add aseptically to the **Growth Factor Mixture (GFM80005)**. Mix well and add asepti-cally to the medium. Mix well before pouring.

Prepared media	
Bottled media:	100 ml: MHT30100, 500 ml: MHT30500
Plated Mueller-Hinton II Agar:	90 mm Petri-dish, 25 ml: MHT50090-01
Plated Mueller-Hinton II Blood Agar, EUCAST:	90 mm Petri-dish, 25 ml: MHT50090-04
Plated Mueller-Hinton II Chocolate Agar:	90 mm Petri-dish, 25 ml: MHT50090-03
Colour of blood free agar:	Yellowish
Colour of blood agar:	Ruby red
Colour of chocolate agar:	Chocolate brown
pH (25 °C):	7,2 - 7,4

**Direction:** If necessary, blood may be added to the melted bottled media according to the direction of the dehydrated media. Dispense aseptically into sterile Petri-dishes. Media in Petri-dishes are ready to use.

## FORMULA in g/l

Acid hydrolisate of casein	17,5
Beef extract	2,0
Starch soluble	1,5
Agar	17,0

**Note:** The typical formula can be adjusted to obtain optimal performance.

**Storage conditions:** Store the dehydrated media tightly closed in a dry place at room temperature. Store the bottled media protected from light at room temperature. Store the plated media protected from light at 2-8 °C. Use before the expiry date on the label.



 $\epsilon$ 

## **Quality control:**

Conditions, Mueller-Hi	nton II Agar:				
Incubation temperature:	J	35 °C	Incubation time:		16 h
Test strains			Growth	Zone diame	eter
Escherichia coli	ATCC25922		Good		
Ampicillin	10 μg			15 – 22 m	ım
Gentamicin	10 μg			19 - 26 m	m
Tigecycline	15 μg			20 – 27 m	ım
TrimethSulfam. 1,25/23	3,75 μg			23 – 29 m	ım
Enterococcus faecalis	ATCC29212		Good		
TrimethSulfam. 1,25/23	TrimethSulfam. 1,25/23,75 μg			26 - 34 m	ım
Pseudomonas aerugino	sa ATCC27853		Good		
Gentamicin	10 μg			17 - 23 m	ım
Tobramycin	10 μg			20 – 26 m	m
Staphylococcus aureus	ATCC29213		Good		
Oxacillin	1μg			19 – 25 m	ım
Gentamicin	10 μg			19 - 25 m	m
Cefoxitin	30 μg			24 - 30 m	ım
TrimethSulfam. 1,25/23	3,75 μg			26 – 32 m	ım

Conditions, Mueller-Hinton II Blood Agar, EUCAST:					
Incubation temperature:	35 °C	Incubation time:	16	h	
Test strains		Growth	Zone diameter		
Streptococcus pneumoniae ATCC49619		Good			
Oxacillin 1 µg			8 – 14 mm		
TrimethSulfam. 1,25/23,75 μg			18 - 26 mm		

Conditions, Mueller-Hinton II Chocolate Agar:					
Incubation temperat	ure:	35 °C	Incubation time:		16 h
Test strains			Growth	Zone diar	neter
Haemophilus influe	nzae ATCC49766		Good		
Ciprofloxacin	5 μg			32 - 40	mm
Cefotaxime	5 μg			29 - 37	mm

**References:** Mueller and Hinton (1941) Proc. Soc. Exp. Biol. Med. 48: 330. www.eucast.org

In vitro diagnostic - for professional use only!

according to Regulation (EC) No. 2015/830

# 1. IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING

1.1 Product identifiers

Product name (dehydrated media):

Product name (bottled media):

Product name (plated media):

Mueller-Hinton II Agar Mueller-Hinton II Agar Mueller-Hinton II Blood Agar,

EUCAST MHT

**Code number:** 

Registration number: - CAS number: -

1.2 Relevant identified uses of the substance or mixture and uses advised against

Recommended use: Microbiological diagnostics
Uses advised against: No information available

1.3 Details of the supplier of the safety data sheet

Manufacturer/distributor: BIOLAB Inc.

H-1141 Budapest, Öv u. 43.

Hungary

Telephone: +36-1-221-9614
Fax: +36-1-364-2006
E-mail: export@biolab.hu

1.4 Emergency telephone number

Emergency telephone: Please contact the regional Authority in your country.

#### 2. HAZARDS IDENTIFICATION

#### 2.1 Classification of substance or mixture

Not a hazardous substance or mixture according to Regulation (EC) No 1272/2008. Not a hazardous substance or mixture according to EC-directives 67/548/EEC or 1999/45/EC.

#### 2.2 Label elements

Labelling in accordance with Regulation (EU) No. 1272/2008:

Signal word: -

Hazard Statements (H phrases): -

Precautionary Statements (P phrases): -

Pictograms: -

#### 2.3 Other hazards

No information available

#### 3. COMPOSITION/INFORMATION ON INGREDIENTS

#### 3.1 Substances

#### 3.2 Mixtures

It is not necessary to publish the components.

according to Regulation (EC) No. 2015/830

#### 4. FIRST AID MEASURES

#### 4.1 Description of first aid measures

General advice: Consult a physician. Show the label and this safety data sheet to

the doctor in attendance. Never administer anything orally to

persons who are unconscious.

After inhalation: Remove the victim from exposure and move into open air.

Consult a physician.

After skin contact: Remove contaminated clothing. Flush the skin with water, then

wash thoroughly with soap and water.

After eyes contact: Rinse out eyes with plenty of clean and cold water while

pulling eyelids up, and seek medical assistance.

After ingestion: Rinse the mouth with water and seek immediate medical

attention. Never administer anything orally to persons who are

unconscious. Never induce vomiting.

#### 4.2 Most important symptoms and effect, both acute and delayed

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

#### 4.3 Indication of any immediate medical attention and special treatment

needed No information available

#### 5. FIRE FIGHTING MEASURES

#### 5.1 Suitable extinguishing media

Use water spray, dry chemical or carbon dioxide.

## 5.2 Special hazards arising from the substance or mixture

In case of fire: Carbon-oxides, nitrogen oxides.

Thermal decomposition can lead to release of irritating gases and vapore.

#### **5.3 Advice for fire-fighters**

Wear self contained breathing equipment for firefighting.

#### 6. ACCIDENTAL RELEASE MEASURES

#### 6.1 Personal precautions, protective equipment and emergency procedures

Avoid generation of dusts; do not inhale dusts. Provide appropriate exhaust ventilation. Use personal protective clothing.

### **6.2 Environmental precautions**

Prevent the contamination of drains, surface or subterranean waters, and the ground.

#### 6.3 Methods and material for containment and cleaning up

Take up dry. Keep in suitable, closed containers for disposal. Clean up affected area.

#### **6.4 Reference to other sections**

Refer to protective measures listed in Section 8 and 13.

#### 7. HANDLING AND STORAGE

#### 7.1 Precautions for safe handling

No special requirement.

Avoid direct contact with the material, spillage, ingestion, eye and skin contact, inhalation. No eating and smoking during working. In the application area, smoking, eating and drinking must

according to Regulation (EC) No. 2015/830

be prohibited. Bathing facilities with hot water, emergency douche and eye irrigator have to be ensured. Wash hands after working with the product.. Remove contaminated clothing.

#### 7.2 Conditions for safe storage, including any incompatibilities

Store the dehydrated media tightly closed in a dry place at room temperature. Store the bottled media protected from light at room temperature. Store the plated media protected from light at  $2-8\,^{\circ}\text{C}$ .

#### 7.3 Specific end use(s)

Use in laboratories.

#### 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

#### **8.1 Control parameters**

The product does not contain substances with Professional Exposure Environmental Limit Values.

**Exposure limits** 

Allowable workplace concentration ( AK ) ( $mg/m_3$ ): No information available Workplace peak concentration: (CK)  $mg/m_3$ ): No information available Biological exposure: No information available

### **8.2 Exposure controls**

General instructions

Avoid direct contact with the material, spillage, ingestion, eye and skin contact, inhalation. No eating and smoking during working. Bathing facilities with hot water, emergency douche and eye irrigator have to be ensured. Wash hands after working with the supplement. Remove contaminated clothing.

Personal protection: Wear protection clothing, remove the contaminated chlothing.

Wash hands after working with the substance.

Personal protection equipment: respirator required (type P1) when dusts are generated

Eye protection: safety glasses are required

Hand protection: use disposable gloves (nitrile caoutchouc 0,11 mm,

breakthrough time 0,11 mm, breakthrough time 480

min.

and wash hands after working with this substance

#### 9. PHYSICAL AND CHEMICAL PROPERTIES

#### 9.1 Information on basic physical and chemical properties

Appearance (dehydrated medium): Homogeneous hygroscopic powder

Colour (dehydrated medium): Yellowish

Odor (dehydrated medium): No information available

pH before autoclaving (25 °C): 7,3 approx.

Appearance (ready to use medium): Gel

Colour (ready to use medium): Yellowish Colour (plated blood agar): ruby red

Odor (ready to use medium): No information available

pH (ready to use medium, 25 °C): 7.2 - 7.4 Water solubility: > 38 g/l

Odor threshold:

Melting point/range:

Flash point:

Evaporation rate:

Flammability (solid, gas):

No information available

Vapor pressure: No information available

according to Regulation (EC) No. 2015/830

Vapor density:

Relative density:

No information available
Partition coefficient (n-octanol/water): No information available
Autoignition temperature:

No information available
Decomposition temperature:

No information available
Viscosity:

No information available
Explosive properties:

No information available
No information available
No information available

**9.2 Other information**No information available

#### 10. STABILITY AND REACTIVITY

#### 10.1 Reactivity

No information available

#### **10.2 Chemical stability**

Stable under the recommended handling and storage conditions (see section7).

#### 10.3 Possibility of hazardous reactions

No information available

#### 10.4 Conditions to avoid

Strong heating, in case of powdered media: risk of powder explosion.

## 10.5 Incompatible materials

No information available

#### 10.6 Hazardous decomposition products

See details in section 5.

#### 11. TOXICOLOGICAL INFORMATION

#### 11.1 Information on toxicological effects

Acute toxicity:

Skin irritation/corrosion:

Serious eye damage/irritation:

Respiratory or skin sensitization:

Germ cell mutagenicity:

No information available
No information available
No information available

Carcinogenicity: No component of this product present at levels greater

than or equal to 0.1% is identified as probable, possible or

confirmed human carcinogen by IARC.

Reproductive toxicity:

STOT – single exposure:

STOT – repeated exposure:

Aspiration hazard:

No information available
No information available
No information available

#### 12. ECOLOGICAL INFORMATIONS

#### 12.1 Toxicity

No information available

#### 12.2 Persistence and degradability

No information available

according to Regulation (EC) No. 2015/830

#### 12.3 Bioaccumulative potential

No information available

### 12.4 Mobility in soil

No information available

#### 12.5 Results of PBT and vPvB assessment

No information available

#### 12.6. Other adverse effects

Do not dispose of the residue into water, sewage or soil.

#### 13. DISPOSAL CONSIDERATIONS

#### 13.1 Waste treatment methods Product:

Residue is classified as hazardous.

**Contaminated packaging:** Handle contaminated empty packaging in the same way as the substance itself.

**Waste treatment regulations**: Dispose of in accordance with the European Directives on waste and hazardous waste. Dispose of in accordance with local regulations.

#### 14. TRANSPORT INFORMATION

#### 14.<u>1</u> UN-number

ADR/RID: - IMDG:- IATA:-

## 14.2 UN proper shipping name

ADR/RID: Not dangerous goods IMDG: Not dangerous goods IATA: Not dangerous goods

#### 14.3 Transport hazard class(es)

ADR/RID: no IMDG: no IATA: no

14.4 Packaging

ADR/RID: no IMDG: no IATA: no

14.5 Environmental Hazards

ADR/RID: no IMDG: no IATA: no

#### 14.6 Special precautions for user

No special precautions required

#### 14.7 Transport in bulk according to Annex II of MARPOL and the IBC code

Not applicable, packaged goods

#### 15. REGULATORY INFORMATION

# 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

Regulation (EC) No 1907/2006 (REACH) and 1272/2008 EK (CLP-GHS)

#### 15.2 Chemical safety assessment

There has been no evaluation a chemical safety assessment of the product.

according to Regulation (EC) No. 2015/830

During recent review the Safety Data Sheet was revised completely.

The information in this Safety Data Sheet is based on current knowledge and on current EC and national laws, as far as the working conditions of the users is beyond our knowledge and control. The product must not be used for purposes other than those that are specified without first having written instructions on how to handle. It is always the responsibility of the user to take the appropriate measures in order to comply with the requirements established by current legislation. The information contained in this

Safety Sheet only states a description of the safety requirements for the preparation, and it must not be considered as a guarantee of its properties.

This document is produced electronically and is valid without signature.



H-1141 Budapest Öv u. 43. Phone: +36 221 9614 Fax: +36 364 2006 e-mail: biolab@biolab.hu, export@biolab.hu

web: www.biolab.hu

2018.05.02. 13:40:43

## **QUALITY CONTROL CERTIFICATE**

# MINŐSÉGI BIZONYÍTVÁNY

Mueller-Hinton II Agar  Mueller-Hinton II agar	
Catalogue number: Katalógus szám:	MHT20500
Packaging (in case of KITs only): Csomagolás (csak KIT-ek esetén):	-
Lot: Lot:	MHT010318125
Manufacture date: Gyártási dátum:	2018. március 01. / 01-03-2018
Expiry: Lejárat:	2021. március / 03-2021

PHYSICAL CO	NTROL					
Fizikai ellenőrze	és					
Dehydrated Media	Criteria				Result	
Portáptalaj	Előírás				Eredmény	
Colour: Szín:	<b>Yellowish</b> Sárgás				Conform Megfelel	
Appearance Megjelenés:		Homogeneous hygroscopic powder Homogén nedvszívó por			Conform Megfelel	
Prepared Media	Sterilization:	121 °C 15 min	Supplement			
Késztáptalaj	Sterilizálás:		Kiegészítő:			
Directions Bemérés: 38 g/l						
	<b>Colour</b> Szín	<b>Clarity</b> Tisztaság	<b>Deposit</b> Üledék		<b>Gel strength</b> Gélerősség	Result Eredmény
Bottled Media	Yellowish	Transparent	None		Suitable for good inoculation	
Késztáptalaj	Sárgás	Transzparens	Mentes		Jól inokulálható	Megf
Plated Media Tubed Media	:	:	:		·	
Csöves táptalaj						4
Lemeztáptalaj						

H-1141 Budapest Öv u. 43.

Phone: +36 221 9614 Fax: +36 364 2006 e-mail: biolab@biolab.hu, biolab-ferol@biolab.hu

web: www.biolab.hu



CHEMICAL CONTROL				
Kémiai ellenőrzés				
Prepared Media	Criteria	F	Result	
Késztáptalaj:	Előírás	E	Eredmény	
pH before autoclaving (25 °C)	7,3 +/- 0,2	7	7,4	
pH autoklávozás előtt (25 ˚C -on)	, ,		<u> </u>	

MICROBIOLOGICAL CON  Mikrobiológiai ellenőrzés			
Control strains Teszt törzsek	Incubation temperature: Inkubációs hőmérséklet:	37 °C	Incubation time: Inkubációs idő:  16 h
Escherichia coli ATCC 25922	Good growth, zone diameters according to EUCAST  Jó növekedés, kioltási zónák EUCAST szerint		Conform Megfelel
Enterococcus faecalis ATCC 29212	Good growth, zone diameters according Jó növekedés, kioltási zónák EUCAST sz		<b>Conform</b> Megfelel
Pseudomonas aeruginosa ATCC 27853	grand, and a second of the sec		Conform Megfelel
Staphylococcus aureus ATCC 29213	Good growth, zone diameters according Jó növekedés, kioltási zónák EUCAST sz	_	<b>Conform</b> Megfelel

H-1141 Budapest Öv u. 43.

Phone: +36 221 9614 Fax: +36 364 2006

e-mail: biolab@biolab.hu, biolab-ferol@biolab.hu





CHEMICAL CONTROL		
Kémiai ellenőrzés		
Prepared Media	Criteria	Result
Késztáptalaj:	Előírás	Eredmény
pH before autoclaving (25 °C)	7,3 +/- 0,2	7,4
pH autoklávozás előtt (25 °C -on)	, ,	,

MICROBIOLOGICA			
Mikrobiológiai ellenő	rzés		
Control strains	Incubation temperature:	0	Incubation time:
Teszt törzsek	Inkubációs hőmérséklet:	37 C	Inkubációs idő:

#### Hereby we declare:

This product has been tested by Quality Control Laboratory and conforms to the specification contained in the relevant catalogue or to the specification agreed with the customer. This product was manufactured by Biolab Inc.

Under our sole responsibility that the above mentioned devices meet the applicable provisions of the Directive 98/79/EC on "In Vitro Diagnostic Medical Devices". All the supporting documents, as required by Annex III on the 98/79/EC Directive, in order to prove conformity to the Essential Requirements as listed in Annex I, are retained under the premises of the Manufacturer.

#### Kijelentjük, hogy

Kizárólagos felelősséget vállalunk, hogy a fenti termék megfelel az EU 98/79 - es In vitro diagnosztikai orvostechnikai eszközökről szóló direktívájában foglalt valamennyi rendelkezésének. Az ezt alátámasztó dokumentáció, amint azt a direktíva III. sz. függeléke megköveteli az I. sz. függelékben felsorolt elsődleges követelményeknek való megfelelés bizonyítása céljából, a gyártónál elérhető.

A terméket a minőségellenőrző laboratórium a fent előírt követelményeknek megfelelőeknek találta. A termék gyártója a BIOLAB Zrt.

This certificate was issued electronically and is valid without signature.

Ez egy elektronikusan készült bizonylat, mely aláírást nem igényel.

Printed by SAP Business Portal

SAP Partner: Mosaic BS Kft

Page 3 of